

Nucleus Smart App Operating System Requirements



Made for
 iPhone | iPad | iPod

- iOS 10.0 and above
- Pair within Settings
- Control with App on iPod or iPhone*

How to check iOS version: Settings →
General → About → Version



- Lollipop 5.0 and Bluetooth 4 and above**
- Pair within the App
- Control with App

How to check Android version: Apps →
Settings → About phone (or About device)
→ Android version

*For a list of compatible iOS devices visit: www.cochlear.com/nucleus/compatibility
**For a list of verified Android devices visit <http://www.nucleussmartapp.com/android>.

Android Phones verified for Nucleus Smart App



Samsung Galaxy S8



Samsung Galaxy S7



Huawei P8lite



**Lenovo
MOTO G5**



Google Pixel



LG G6



HTC U11

